# Analytical Semi-Micro Balances



High Performance of up to 0.01mg Semi-Micro Weighing Automatic Self-Calibration and Ergonomic Design







# Superb "New Generation

Speedy & Simple Operation to Maximize Performanace

# **Easy Access Door**

A unique system allows the user to conveniently open or close the left or right side of the weighing chamber using a control lever.

This system allows quick and efficient access to the weighing chamber.





#### **One Touch Calibration**

Allows quick calibration with one push of the CAL key.

### **Automatic Adjustable Environment Setting**

Adapts automatically to the most suitable Environment Response Setting to minimize adverse external effects. "FAST", "MID" or "SLOW" is displayed.



### **Automatic Self Calibration**

Performs self calibration automatically when a change of ambient temperature exceeds the preset temperature range.



#### 13 Multiple Weighing Units

g, mg, pcs, %, oz, ozt, ct, mom, dwt, GN, tl, t(tola), m (messghal)

#### Two Layer Housing

Durable construction minimizes negative effects of ambient temperature changes.

## Splash Proof Keyboard and Display

Protected from dust and spills.

#### **Data Memory Function**

Stores up to 200 sets of weighing data. Data can be output using numerical codes.



#### **GLP/GMP Compliant**

Allows GLP or LIMS balance management by outputting the Balance ID number and weight used to calibrate the balance. The data can be output to A&D's AD-8121B printer or a computer, indicating date, time, Balance ID number and calibration weight.





















# n" Analysis!



The GR Series has a 30% smaller footprint compared with conventional balances.

Its tall weighing chamber can accommodate large beakers.



#### Standard Underhook

Ideal for density measurement and weighing magnetic substances

#### Standard RS-232C

Enables bi-directional communication with a PC or printer.

## **Quick Reference Card**

A fast and convenient operation guide.

#### **WinCT-Window Communication Tools**

WinCT is Windows Communication Tools Software, providing easy transmission of weighing data from A&D Balances to a PC using RS-232C interface.


Data can be transferred in a manageable format limited only by a PC's capability.

#### WinCT consists of three applications:

**RsKey** allows weighing data to be directly pasted into a Windows application such as EXCEL or WORD, but this application cannot use "commands" to control the balance.

**RsCom** can transmit weighing data to a PC using a text file format and can transmit "commands" to control the balance. "Commands" include Calibrate, On/Off, Print, Query weight data, Re-zero, Range, Sample, Mode and Zero.

**RsWeight** can retrieve weighing data from a balance and display it in graph form on a PC screen in realtime. Maximum, minimum, average, standard deviation and coefficient of variation values of data can be calculated and displayed.



RsCom, RsKey and RsWeight offer user selectable time, date, sequential #, interval transmission and GLP/GMP data. A readme file contains all setup, operation & troubleshooting guides.

#### **Specifications**

			GR-202	GR-300	GR-200	GR-120
Gram	(g)		210 x 0.0001 / 42 x 0.00001	310 x 0.0001	210 x 0.0001	120 x 0.0001
	(mg)		210000 x 0.1 / 42000 x 0.01	310000 x 0.1	210000 x 0.1	120000 x 0.1
Decimal Ounce (oz)			7.4 x 0.00001 / 1.4 x 0.000001	10 x 0.00001	7.4 x 0.00001	4.2 x 0.00001
Troy Ounce (ozt)			6.7 x 0.00001 / 1.3 x 0.000001	9.9 x 0.00001	6.7 x 0.00001	3.8 x 0.00001
Pennyweight (dwt)			135 x 0.0001 / 27 x 0.00001	199 x 0.0001	135 x 0.0001	77 x 0.0001
Carat (ct)			1050 x 0.001 / 210 x 0.0001	1550 x 0.001	1050 x 0.001	600 x 0.001
Momme (mom)			56 x 0.0001 / 11 x 0.00001	82 x 0.0001	56 x 0.0001	32 x 0.0001
Grain Unit (GN)			3240 x 0.002 / 648 x 0.0002	4784 x 0.002	3240 x 0.002	1851 x 0.002
Tola (t)			18 x 0.00001 / 3.6 x 0.000001	26 x 0.00001	18 x 0.00001	10 x 0.00001
Tael (TL)			5.5 x 0.00001 / 1.1 x 0.000001	8.2 x 0.00001	5.5 x 0.00001	3.1 x 0.00001
Percentage Min Div.			0.01%			
Counting Min Weight			0.1mg			
Linearity			±0.0002 / ±0.00003g	±0.0003g	±0.0002g	
Repeatability/Std. Dev.			0.0001 / 0.00002g	0.0002g	0.0001g	
Stabilization Time			3.5 / 8 seconds (typically)	ds (typically) 3.5 seconds (typically)		
Sensitivity Drift			±2ppm / °C (10°C~30°C / 50°F~86°F)			
Display Refresh			5 times per second / 10 times per second			
Pan Size			Ø85mm / 3.3inches			
Physical I	Dimensions	mm	249 (W) x 330 (D) x 328 (H)			
		inches	9.8 (W) x 12.99 (D) x 12.91 (H)			
Breeze Bre	eak Dimensions	mm	178 (W) x 160 (D) x 233 (H)			
		inches	7.01 (W) x 6.5 (D) x 9.17 (H)			
Admissible Ambient Conditions			5°C~40°C / 41°F~104°F RH less than 85%			
Weight Approximately			6.0kg / 13.2lb			
Power			Approx. 11VA (supplied to the AC adapter)			
Standard Accessories			WinCT in CDROM, Manual, AC Adapter, Reference Card			

Specifications are subject to change for improvement without notice

#### **Accessories**

**AD-1653** Density Determination Kit **AD-8121B** Compact Printer

AD-1651 Vibratory Spoon
AD-1682 Rechargeable Battery Unit

AD-1683 Static Eliminator
AD-8920 Remote Display
AD-8922 Remote Controller


**AD-8526** RS/Ethernet Converter with WinCT-Plus


**AX-K01710** RS-232C Cable (25P-9P)

#### AD-1653 / Optional Density Determination Kit





#### GR Series Dimensions (mm/inches)





#### -----







#### Connect With Us













