


Microtome knife sharpener type MMS

Safe and easy to use automatic
microtome knife sharpener



Microtome knife sharpener type MMS




Safe and easy to use automatic
microtome knife sharpener

SLEE have produced and reconditioned steel and hard-metal (tungsten-carbide) knives for decades in the factory in Mainz. However, it is not always economical to find a suitable factory for this service in your area or return the knives for reconditioning to us. Therefore this low-cost microtome knife sharpener type MMS, is a valuable alternative and addition to our manufacturing programme

- ▶ For microtome knives profile C and D with a length of 100 – 250 mm
- ▶ Variable honing programs by means of: 2 honing compounds
- ▶ 2 honing surfaces (coarse + fine)
- ▶ Variable speed, selectable
- ▶ Time clock up to 120 minutes
- ▶ Height adjustment for better sharpening result
- ▶ High effectiveness by means of glass-plate oscillation and selection of honing
- ▶ Program depending on amount of wear of knife
- ▶ Sight check during honing
- ▶ Safety switch stops unit as soon as plastic cover is opened
- ▶ Universal power input 100 – 240 V, 50/60 Hz
- ▶ Accessories:
 - Knifeholder (standard – for knives with a back < 10 mm)
 - Glass honing plate
 - Bottle honing compound fine, 50 ml
 - Bottle honing compound coarse, 50 ml




Sight check during honing



Profile A

Profile B



Profile C

Profile D

| Weight and Dimensions | Depth x Width x Height |
|--|---------------------------------|
| Dimensions (in mm) | 380 x 370 x 280 |
| Dimensions packed (in mm) | 550 x 650 x 480 |
| Weight net | 17 kg |
| Weight packed | 29 kg |
| Electrical parameters | Art.-No. 10200000 230 V / 50 Hz |
| | Art.-No. 10200001 230 V / 60 Hz |
| | Art.-No. 10200002 115 V / 60 Hz |
| | Art.-No. 10200003 100 V / 50 Hz |
| | Art.-No. 10200004 100 V / 60 Hz |
| Please ask for other units for Histopathology and consumables as well as laminar flow equipment. | |



Connect With Us

