# 3011 / 3012 📼

#### Gentle and quiet

Gentle and uniform rocking motion make this quiet unit stand out in every-day use. Both of these units are maintenance-free, carry the GS safety approval and come with the CE mark.



#### Technical Data

Overall dimensions (W x D x H): 510 x 625 x 168 mm

Moving platform: 450 x 450 mm

Load: max.15 kg

Shaking motion: three-dimensional orbital rocking

#### Timer:

3011 / up to 60 minutes, or continuous operation 3012 / 1 min. - 99:59 hours, or continuous operation

Shaking amplitude: 3 degrees

Shaking frequency: 2 - 50 rpm

Electrical connection: 230 V\* \* Other voltages on request! 3011 / 50...60 Hz, 90 W 3012 / 50 or 60 Hz\*\*, 90 W \*\* Please indicate when ordering!

Net/gross weight: 18/22 kg

Packing volume (cardboard box): 0.19 m  $^{\scriptscriptstyle 3}$ 

 Order No. 3011
Order No. 3012 with RS 232 port



**3011** Analogue Orbital Rocking Motion Shaker with Non-slip Mat 3965

#### Features

- compact, low-wear counterbalanced drive mechanism, ensuring dependability for continuous operation with three-dimensional motion
- housing made of high impact strength polystyrene, off-white varnish. The base plate is made of electrolytically galvanised, powdercoated sheet steel, the shaking platform of anodised aluminium, equipped with four plastic pins for secure attachment of accessory equipment
- an a.c. motor with overload protection drives the unit
- 3011 clearly laid-out control panel for easy operation
- 3012 speed and remaining time indicated by two LC displays

- 3011 electronic speed control, stepless, gentle start-up
- 3012 microprocessor-controlled, adjustable in steps of 0.5 rpm, gentle start-up
- constant speed during continuous operation, independent of load

#### Applications

In AIDS research, for diagnostic tests with test plates, or for growing cell cultures and micro-organisms this unit has become indispensable.

Also well-suited for use in laboratories, incubation rooms and moderating rooms at ambient temperatures between +10 °C and +50 °C.

## Accessories



3011



3012

3014

Shaking tray

Clamps

materials

made of stainless steel,

accommodate clamps for

accessory equipment

Erlenmeyer flasks and other

made of stainless steel for Erlenmeyer flasks, for screwing on to shaking tray 3966, complete with attachment

450 x 450 mm, with holes to



3013



Order No. 3966



Order No. 3983	B for	25	ml	flasks	(7	9*)
Order No. 3984	l for	50	ml	flasks	(4	9*)
Order No. 3985	<b>5</b> for	100	ml	flasks	(3	6*)
Order No. 3986	<b>b</b> for	200	ml	flasks	(2	2*)
Order No. 3987	for	250-300	ml	flasks	(1	6*)
Order No. 3988	<b>B</b> for	500	ml	flasks	(1	2*)
Order No. 3989	for	1000	ml	flasks	(	9*)
Order No. 3990	) for	2000	ml	flasks	(	4*)
Other clamps on request						

\* Max. number of clamps per shaking tray





Order No. 3965

#### Non-slip mat

for the shaking platform, 420 x 420 mm size, for slow moving, e.g. nutrient solutions in Petri dishes

Order No. 3967

Universal mount for secure attachment of different shaking objects between the six rubbercoated bars



Order No. 3968

#### Platform frame

with four levels, with top rims, three levels removable. Frame and trays made of stainless steel. Distance between trays - 60 mm. Grooved rubber mats on trays to ensure secure grip of objects.

### Connect With Us f У 🖸 🚫 💶 🔞 🚥