

## VWR® Combi-Syringes for Eppendorf® Repeater® Pipettor



Designed exclusively for VWR repeating pipettors and Eppendorf Repeater pipettors. Manufactured from virgin polypropylene with a polyethylene inner plunger. Each reservoir tip is graduated to indicate the total volume and the repetitive volumes that it can dispense.

Note: Combi-syringes will not fit the Eppendorf Repeater Plus, the Eppendorf Repeater Pro pipettors, or the Eppendorf Repeater Stream pipettors.

### Description

Accuracy	Pipetting Volumes	Cat. No.	Precision	Volume	Unit
Sterile Combi-Syringes	±2.5%	1 µL, 2 µL, 3 µL, 4 µL, 5 µL	40000-056	±5.0–2.5%	0.05 mL Case of 100
Nonsterile Combi-Syringes	±2.5%	1 µL, 2 µL, 3 µL, 4 µL, 5 µL	40000-058	±5.0–2.5%	0.05 mL Case of 100
Sterile Combi-Syringes	±0.8%	10 µL, 20 µL, 30 µL, 40 µL, 50 µL	40000-012	±0.7–0.4%	0.5 mL Case of 100
Nonsterile Combi-Syringes	±0.8%	10 µL, 20 µL, 30 µL, 40 µL, 50 µL	40000-000	±0.7–0.4%	0.5 mL Case of 100
Sterile Combi-Syringes	±0.8%	25 µL, 50 µL, 75 µL, 100 µL, 125 µL	40000-014	±0.5–0.2%	1.25 mL Case of 100
Nonsterile Combi-Syringes	±0.8%	25 µL, 50 µL, 75 µL, 100 µL, 125 µL	40000-002	±0.5–0.2%	1.25 mL Case of 100
Sterile Combi-	±0.7%	50 µL, 100 µL, 150 µL, 200	40000-016	±0.4–0.2%	2.5 mL Case of 100

Syringes		μL, 250 μL				
Nonsterile Combi-Syringes	±0.7%	50 μL, 100 μL, 150 μL, 200 μL, 250 μL	40000-004	±0.4–0.2%	2.5 mL	Case of 100
Sterile Combi-Syringes	±0.4%	100 μL, 200 μL, 300 μL, 400 μL, 500 μL	40000-018	±0.3–0.2%	5.0 mL	Case of 100
Nonsterile Combi-Syringes	±0.4%	100 μL, 200 μL, 300 μL, 400 μL, 500 μL	40000-006	±0.3–0.2%	5.0 mL	Case of 100
Sterile Combi-Syringes	±0.3%	250 μL, 500 μL, 750 μL, 1000 μL, 1250 μL	40000-020	±0.3–0.1%	12.5 mL	Case of 100
Nonsterile Combi-Syringes	±0.3%	250 μL, 500 μL, 750 μL, 1000 μL, 1250 μL	40000-008	±0.3–0.1%	12.5 mL	Case of 100
Sterile Combi-Syringes	±0.3%	500 μL, 1000 μL, 1500 μL, 2000 μL, 2500 μL	40000-054	±0.2–0.1%	25.0 mL	Case of 100
Nonsterile Combi-Syringes	±0.3%	500 μL, 1000 μL, 1500 μL, 2000 μL, 2500 μL	40000-052	±0.2–0.1%	25.0 mL	Case of 100
Sterile Combi-Syringes	±0.2%	1000 μL, 2000 μL, 3000 μL, 4000 μL, 5000 μL	40000-022	±0.2–0.1%	50.0 mL	Case of 100
Nonsterile Combi-Syringes	±0.2%	1000 μL, 2000 μL, 3000 μL, 4000 μL, 5000 μL	40000-010	±0.2–0.1%	50.0 mL	Case of 100
Combination Set	-	-	40000-050	-	-	Case of 100