



Chromatography Column with Flow Control Connection, Ace Glass Incorporated




This chromatography column has an integral reservoir for connecting flow control adapters such as 89053-324 ([SJ] 35/20) and 89053-326 ([SJ] 50/30). This column has a 2mm bore, 1:5 taper solid PTFE stopcock plug at bottom. It is available with or without a Porosity B (70–100 micron) fritted disc sealed in. Fabricated from borosilicate glass.

Ordering Information: Pinch clamp (89060-258) and flow control adapter (89053-324 for [SJ] 35/20 and 89053-326 fo [SJ] 50/30) must be ordered separately.




35/20 10.5MM POR B COLUMN



VWR Catalog Number	Quantity	Description	[ST] Joint Size	Capacity	Reservoir Capacity	Supplier No.
89053-598		Chromatography Column With 70-100 micron Fritted Disc, Fits Tubing ID 10.5 mm, 250 mm	35/20	22 mL	200 mL	5906-33



89053-600		Chromatography Column With 70-100 micron Fritted Disc, Fits Tubing ID 14.5 mm, 250 mm	35/20	41 mL	250 mL	5906-36
-----------	---	---	-------	-------	--------	---------



89053-602				Chromatography Column With 70-100 micron Fritted Disc, Fits Tubing ID 19.0 mm, 300 mm	35/20	85 mL	250 mL	5906-39
GSA (SBE)								
89053-604				Chromatography Column With 70-100 micron Fritted Disc, Fits Tubing ID 25.0 mm, 300 mm long	35/20	147 mL	250 mL	5906-43
GSA (SBE)								
89053-606				Chromatography Column With 70-100 micron Fritted Disc, Fits Tubing ID 38.0 mm, 300 mm long	50/30	340 mL	250 mL	5906-45
GSA (SBE)								
89053-608				Chromatography Column With 70-100 micron Fritted Disc, Fits Tubing ID 50.0 mm, 300 mm	50/30	589 mL	500 mL	5906-47
GSA (SBE)								
89053-586				Chromatography Column Without Disc, Fits Tubing ID 10.5 mm, 250 mm long	32/20	22 mL	200 mL	5906-06
GSA (SBE)								
89053-588				Chromatography Column Without Disc, Fits Tubing ID 14.5 mm, 250 mm long	35/20	41 mL	250 mL	5906-09
GSA (SBE)								
89053-590				Chromatography Column Without Disc, Fits Tubing ID 19.0 mm, 300 mm	35/20	85 mL	250 mL	5906-11
GSA (SBE)								
89053-592				Chromatography Column Without Disc, Fits Tubing ID 25.0 mm, 300 mm long	35/20	147 mL	250 mL	5906-14
GSA (SBE)								
89053-594				Chromatography Column Without Disc, Fits Tubing ID 38.0 mm, 300 mm long	50/30	340 mL	250 mL	5906-16

Connect With Us

