

## Lab-Tek\* II Chamber Slide\* System, Sterile, NUNC\*



This system features removable, polystyrene chambers that are sealed to a 25x75mm non-fluorescent glass slide using biocompatible acrylic adhesive. The sealant remains with the media chamber after detachment, leaving a hydrophobic border to define cell culture areas. A tabbed cover fits over the media chamber allowing gas exchange while minimizing evaporation. The slides are provided sterile with an RS wash to create a consistent growth surface for adherent cell lines.

Numbers along the border of the chamber and a Superfrost\* writing area on the slide facilitate identification.

Working Volume	No. of Chambers	Well Area	Cat. No.
2.0-4.5 mL	1	8.6 cm <sup>2</sup>	62407-290
1.0-2.0 mL	2	4.0 cm <sup>2</sup>	62407-292
0.5-1.0 mL	4	1.7 cm <sup>2</sup>	62407-294
0.2-0.5 mL	8	0.7 cm <sup>2</sup>	62407-296



Connect With Us

