

Haemometers acc. to Sahli

Haemometers are used for the determination of blood's content of haemoglobin. The Marienfeld Superior haemometer according to Sahli is supplied as complete set consisting of:

- polystyrene support with 2 coloured rods and opal glass plate
- comparator tube
- haemoglobin pipette 20 µl acc. to Sahli
- silicone tubing of approx. 16 cm length
- white mouthpiece
- dropping pipette with rubber teat
- stirring rod
- acid vial
- cleaning brush
- directions for use
- without CE-mark, only for sale and application outside EC

Cat. No.	Multipack	Unit
▶ 32 430 00	60 x 1	1

Comparator tubes

- made of soda lime glass
- with round bottom and beaded rim
- with double coloured scale yellow/red
- without CE-mark, only for sale and application outside EC

Cat. No.	Multipack	Unit
▶ 32 440 00	100 x 10	10

Dropping pipettes

- made of soda lime glass
- with straight tip
- with beaded rim
- without rubber teat
- handy packaging in a carton

Cat. No.	Length	Diameter	Unit
▶ 40 100 20	75 mm	6/7 mm	100
▶ 40 100 70	150 mm	7/8 mm	200

Rubber teats for dropping pipettes

Cat. No.	Suction volume	Colour	Unit
▶ 59 020 02	approx. 0.8 ml	red	100



When pipetting, please note the safety regulations. Pipetting by mouth is dangerous. We disclaim any responsibility.

The scales of our comparator tubes are adjusted to our present haemometer's coloured glass rod. The tubes are not to be used for application with haemometers we produced in the past or for those of other manufacturers.

